MEDICO
Italian Design for Life

ITALIANO

FRANÇAIS

ENGLISH

ESPAÑOL

ČESKÝ

SLOVENSKÝ

Vivere con il pacemaker
Vivre avec le stimulateur cardiaque
Living with a pacemaker
Vivir con el marcapasos
Život s kardiostimulátorem
Život s kardiostimulátorom

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The heart

The heart is a muscular organ which continuously supplies the blood to all parts of the body thanks to its cyclical and constant pump activity.

It is composed of four cavities: two atria and two ventricles. The right atrium receives the venous blood from the whole body and pumps it into the right ventricle. Here, through the pulmonary vein the blood is pumped to the lungs where it is recharged with oxygen. The oxygen-rich blood returns to the left atrium, is conveyed to the left ventricle and through the arteries is then distributed to all the tissues.

The cardiac muscular contractions are due to an initial electrical stimulus which is produced by the heart itself.

The sum of all the synchronous contractions of the cardiac muscle fibres is such as to generate a push to the blood contained into the heart chambers. That way the blood is spread across the circulatory system and is therefore able to reach every remote part of the body. Thus, the heart plays a key role. At rest it beats from 60 to 80 times per minute, moving during this period approximately 5-6 litres of blood. The rate of those beats, i.e. the number of contractions per minute, could change according to a body request of a major blood amount, as during physical activity or emotional stress (scares and strong emotions).

Particular cardiac diseases and some illnesses of the heart conduction system could imply some cardiac rate alterations.

The pacemaker

The pacemaker is a device that constantly monitors the patient's cardiac rate and generates electrical impulses when required, i.e. when there is a rhythm slowdown, in order to grant an optimal heart rate. The electrical impulses produced by the pacemaker are conducted to the heart through special leads called electrocatheters. This device could also store an enormous amount of information on the heart activity, which is an useful instrument for the physician for an optimal therapy definition.

The pacemaker does not affect everyday activities and hobbies, does not influence work activity (except for heavy labour and works on environments where strong electromagnetic fields are present) and rather constitutes an additional security and a possibility to improve one's physical condition.

The implantation procedure

The pacemaker application involves surgery, and consequently requires a sterile environment, specific radiological equipment and monitoring systems. The operation is generally performed under local anaesthetic and is not time-consuming (from 30 to 60 minutes, depending on the implantation kind).

The lead is inserted into a vein through an incision a few centimetres long, is then driven to the heart and positioned in such a way as to obtain contact with the inner wall (into the right atrium or the right ventricle).

As the leads good position is confirmed, they are connected to the pacemaker. The pacemaker is then positioned inside the body, into a pocket created by the physician under the skin, where the incision was made.

After implantation

During the postoperative period, it is important to follow all the physician's indications.

For the first next days the patient must be careful to limit his own actions and movements of the shoulder corresponding to the pacemaker implant.

Moreover, for a few days from the operation, the patient may feel a slight pain where the incision was made.

The stitches used to close the incision are removed in the hospital after 8-10

days. In this occasion it is also checked the healing state of the scar.

Checking the pacemaker

Medical check-ups are defined by the physician and it's really important to respect them always. During those follow-ups the physician and the hospital personnel perform a variety of instrumental investigations in order to check the goodness of the device functioning and of the applied therapy.

The follow-ups are regularly performed and consist on the pacemaker interrogation with a specific device called programmer. The programmer allows reading all the information stored in the pacemaker and permits the physician to execute specific tests.

If necessary, the hospital personnel also have the possibility to edit the device settings in order to obtain the best programming for every patient.

So it is really important that the patient respects the pacemaker follow-ups.

The pacemaker Remote Monitoring

Remote monitoring is realized with a device that can interact with some pacemaker models produced by MEDICO S.p.A..

This device can interrogate a pacemaker, read the parameters and information stored in it and automatically create a report containing the programming parameters of the stimulator. It is also possible to send this report to the physician.

For detailed information about the best available options for your clinical condition treatment, please contact your physician.

Pacemaker replacement

During the periodic pacemaker follow-ups, the physician checks the charge level of the battery.

When the battery is starting to run out, in other words when the pacemaker is entering the so-called rundown phase (which can last for months), the physician programs its replacement.

The operation is a short surgical procedure which is performed in day hospital (the patient returns home the same day of the replacement).

The procedure requires a little incision on the original implant site. The old device is disconnected from the leads and is replaced with a new pacemaker. The leads are normally left in place after checking their integrity and correct functioning.

The new pacemaker is then connected to the leads and is inserted in the same pocket of the older stimulator.

After replacement it is important that the patient observe the same precautions he was given at the time of his previous pacemaker implant.

Living with a pacemaker

After the implantation and its recovery, the pacemaker user can get back to his activities as and better than before, especially on those activities that had been compromised by the cardiac dysfunction.

Physical exercise. The pacemaker does not imply any limitations to physical activities.

Physical exercise must be proportionate to the patient's general condition and to his personal resistance limits and exertion adaptability. It is important to avoid bumps, hits and excessive solar exposure on the pacemaker site, during both physical and usual activities. Some normal activities that a pacemaker user can do with any problem are for example: gardening, bicycling, climbing stairs, doing housework, swimming, etc. Sexual activity is not contraindicated. More demanding physical activities must be discussed with the physician. For any doubt or question, please contact your physician.

Diet. No particular restrictions or dietary concerns are required. It is however a good habit to avoid heavy meals and to assume a moderate quantity of coffee and alcohols. Particular cases should either way follow the diet recommended by the physician.

Drugs. The pacemaker does not replace drugs acting on the cardiocirculatory system, but along with them it contributes to the patient's clinical therapy. So, it is important to take all drugs recommended by the physician. The pacemaker is not a contraindication to all other types of drugs.

Journeys. The pacemaker is not a limitation to journeys by car, motorcycle, train, ship and plane, if the patient's clinical condition is good enough. It is a good habit to carry the pacemaker User's Identity Card during journeys, together with the physician's contact information, name and dosage of drugs.

Electromagnetic anti-theft systems

(e.g. on airports, tribunals or banks).

Inform security personnel of the implanted pacemaker. The metal parts of the device may indeed turn on the metal detector alarm.

The pacemaker is protected from common electrical interferences and it could hardly be affected by the interferences induced by anti-theft systems of airports, tribunals, banks, museums, etc. Avoid however standing or leaning against a metal detector. That way it is reduced the possibility of a temporary interference occurrence.

Ask if hand search is possible. If hand probes are used, inform security personnel that the detection device shouldn't be kept over the pacemaker, but that it must be passed quickly over it.

Electrical Equipment

The pacemaker is protected from common electrical interferences and it has been designed to be immune from electrical disturbances that could be emitted from common electrical equipment. Nevertheless some of those devices can emit prolonged interference signals which can temporarily affect the pacemaker behaviour.

Possible sources of interference include: radio transmitters and mobile phones, anti-theft systems, electric current conductors, high voltage systems, and equipment for clinical use normally employed for diagnostic, therapeutic or surgical purposes.

If the patient suspects electrical interference near an electrical device and feels dizzy or has an unexpected increase or an irregularity of his cardiac rate, then he must move away from it.

The following lists may not be complete. For any doubt or question, please contact your physician.

Safe. The following items can be normally used, as long as in good condition and well-functioning: kitchen appliances; household appliances; bathroom appliances (electric toothbrushes, hair dryers); electric heater, portable stoves, electric blankets; AM and FM radios, televisions, game consoles, fax and copy machines, CD and DVD players,

VCRs and video recorders, remote controls.

Pay attention and avoid positioning the following items directly over the pacemaker: headphones, earphones, stereo amplifiers; electric razors and shavers; cordless household telephones.

Pay attention also at:

Cellular Phones. They can be safely used as long as a 15 cm distance is always kept between the phone and the pacemaker. For security reasons it is recommended to make and receive phone calls using the ear at the opposite side to the pacemaker. Also avoid putting the phone near the pacemaker, for example in the shirt or breast pocket over it.

PCs, Tablets and palm PCs, Wi-Fi and Bluetooth devices. The risk involved by those devices varies from case to case. Those devices should be however kept at a distance of at least 15 cm from the stimulator. That way it is assured a wide safety margin.

Security systems in shops. Walk steadily through the anti-theft detectors and avoid lingering in close proximity.

Alternators and motors of cars and any other internal combustion engine. Avoid leaning over, especially if running, and keep a 30 cm distance from the stimulator.

Drills and wire-powered tools,

chainsaws, table saws.

When used, a 30 cm distance must always be kept between those tools and the pacemaker.

Sunlamps and tanning beds.

Avoid overheating of the area where the pacemaker is set.

Possible sources of magneto-static interference (magnets,

electromagnets, magnet poles of electric motors, magnetic tags and signs, etc.). A 15 cm distance is always recommended for everyday objects. However the effect of the magnetic field immediately disappears as the magnet is removed or distanced.

Firearms. If the pacemaker is

set in the same shoulder where the rifle stock is placed, then the use of it is strongly discouraged, as the kickback (recoil) can irreversibly damage the pacemaker and the leads.

Avoid if possible the following items: clothes with magnetic buttons; magnetic mattresses, mats, chairs and cushions; body fat scales and monitors; equipment that produce strong vibrations (jackhammers, pneumatic drills), arc welders and resistance welders; high-power electrical transformers and high-voltage industrial equipment; radio transmitters and radar stations; big electric motors (avoid leaning over, especially if functioning); CB radio antennas and high power radio-telephones, especially the ones with TETRA system; faulty or broken electrical equipment,

with possible short circuits.

Medical and Dental Procedures

It is important that the patient informs of his condition and of the presence of a pacemaker the general practitioner, the dentist, the physiotherapist, the surgeon and the emergency and hospital personnel, every time it could be useful. The hospital personnel must be advised every time the patient has to undergo a medical procedure.

For more information please visit the patient section on:

www.medicoweb.com

This guide to pacemaker users is intended for purely informational purposes, not as a substitute to medical advice nor can it be considered an appropriate instrument of medical consultancy.

For any doubt, question or request of information, please consult your physician.

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